

Alsto LABORATORY FLOTATION TEST APPARATUS

The **Alsto** Laboratory Flotation Test Apparatus is used primarily for metallurgical test work to establish reagent requirements in production equipment.

Features include:

- Variable digital shaft RPM display.
- Air pressure and flow control
- Digital timing.
- Password setting
- 50 & 60 Hz
- Single phase 0.75kW

Options:

Acrylic cells:	0.5; 1.75, 3.5, 5.0 & 7.0 litre
Stainless steel cells:	0.5; 1.75, 3.5, 5.0 & 7.0 litre
PE impellers:	0.5; 1.75, 3.5, 5.0 & 7.0 litre

Can be used as:

- Agitator
- Scrubber
- Attritioner
- Blender
- Mixer



ALSTO

Phone: +61 (08) 9526 2523 | Email: sales@alsto.com.au
Fax: +61 (08) 9526 2522 | ABN: 41 047 362 761

HS Code 847410 & Associated spares 847490